



*Python-based **H**ierarchical **EN**vironment for **I**ntegrated **X**tallography*

## **Tutorial: Fit Biomolecules into Cryo-EM Maps using MD Simulation (commandline)**

### **Contents**

- [Tutorial: Fit Biomolecules into Cryo-EM Maps using MD Simulation \(commandline\)](#)
- [Overview](#)
- [Input files](#)
- [Initial Model](#)
- [Target Map](#)
- [Running the program](#)
- [Output](#)
- [Examples](#)
- [Limitation](#)
- [Reference](#)
- [Options](#)
- [List of most useful options](#)

## **Overview**

This tutorial will show you how to fit biomolecule atomic structures into cryo-EM maps using molecular dynamics simulation with PHENIX commandline

For GUI execution, please see the [cryo\\_fit\\_gui\\_tutorial](#)

Theoretical explanation of cryo\_fit is [here](#)

For installation of cryo\_fit, please see the [installation notes for cryo\\_fit](#)

## **Input files**

<initial\_model> and <target\_map>

## **Initial Model**

Available format: .cif and .pdb

The initial model is a guide or template structure (CIF/mmCIF/pdb) that is close to a target cryo EM map structurally.

You can use either [map\\_to\\_model](#) or UCSF chimera (Tools -> Volume Data -> Fit in Map)) to prepare the initial model.

## **Target Map**

Available format: .ccp4 and .map (MRC style in binary file) and .sit (Situs style in text file)

## **Running the program**

```
% phenix.cryo_fit <initial_model> <target_map>
```

example command line:

```
% phenix.cryo_fit tRNA.pdb tRNA.map
```


## Output

A final cryo\_fitted structure: steps/8\_cryo\_fit/cryo\_fitted.pdb (and cryo\_fitted.gro for vmd visualization)

Correlation coefficients (CC) between cryo\_fitted structures and cryo-EM maps: steps/8\_cryo\_fit/cc\_record

## Examples

cryo\_fit example with tRNA



Before cryo\_fit  
(cc = 0.53)

After cryo\_fit  
(cc = 0.69)

## Limitation

gromacs4.5.5 seems to not handle H2O (water) heteroatom. cryo\_fit will remove water molecules (if any) from the input .cif/.pdb and fit to cryo electron microscopy map.

cryo\_fit doesn't handle non-canonical "residue"s such as 7C4, BMA, GDP, ILX, NAG, SEP, TRX. The cryo\_fit will simply erase those residues.

## Reference

S. Kirmizialtin, J. Loerke, E. Behrmann, C. MT. Spahn, K. Y Sanbonmatsu, Using Molecular Simulation to Model High-Resolution Cryo-EM Reconstructions, Methods Enzymol., 558, 2015, 497-514

## Options

All options will be used as default if unspecified. Gromacs expert users are welcome to customize those options if they wish.

## List of most useful options

Option	Default value	Description of inputs and uses

emweight_multiply_by	8	Multiply by this number to the number of atoms for weight for cryo-EM map bias. For example, emweight = (number of atoms in gro file) x (emweight_multiply_by which is 8) The higher the weight, the stronger bias toward EM map rather than MD force field and stereochemistry preserving constraints. If user's map has a better resolution, higher value of emweight_multiply_by is recommended since map has much information. If user's map has have a worse resolution, lower value of emweight_multiply_by is recommended for more likely geometry. If CC (correlation coefficient) needs to be improved faster, higher number of emweight_multiply_by is recommended.
number_of_cores_to_use	max cores	Specify number of cores for minimization and cryo_fit. If it is not specified, or max is chosen, the cryo_fit will try to use most cores automatically (up to 16)
number_of_steps_for_cryo_fit	None	This is the initial number of steps for cryo_fit. Eventually, cryo_fit will increase it depending on molecule size and cc trend. For tutorial files, this will be 70,000
number_of_steps_for_minimization	None	Specify number of steps for minimization. If this is left blank, cryo_fit will estimate it depending on molecule size.number of steps for cryo_fit. Enough minimization will prevent "blow-up" during MD simulation later.