



How to generate and record movie by cryo_fit1?

Contents

- [How to generate and record movie by cryo_fit1?](#)
- [Pymol \(confirmed with version 2.2.2\)](#)
- [UCSF Chimera \(confirmed with version 1.13\)](#)
- [UCSF ChimeraX \(confirmed with version 0.7 and 0.8 \(2018-08-14\)\)](#)
- [VMD \(confirmed with version 1.9.3\)](#)
- [Note](#)

Pymol (confirmed with version 2.2.2)

[command] cd output/trajectory

[command] load trajectory.gro, mytraj

[command] load_traj traj.xtc, mytraj

Click salmon color play button (triangle) in right bottom to play

(reference) https://pymolwiki.org/index.php/Load_Traj

UCSF Chimera (confirmed with version 1.13)

[menu] Tools -> MD/Ensemble analysis -> MD Movie -> Trajectory format = GROMACS .tpr = (output/trajectory/for_cryo_fit.tpr), .xtc = (output/trajectory/traj.xtc) -> OK

[menu] Tools -> Depiction -> Rainbow -> OK

[menu] File -> Open -> .../user.mrc (set opacity into 50%)

[menu] Presets -> Publication 1

(in MD Movie) click black play button to play

(in MD Movie) File -> Record Movie, to record as .mp4 file

UCSF ChimeraX (confirmed with version 0.7 and 0.8 (2018-08-14))

[menu] File -> Open -> .../output/cryo_fitted_chain_recovered_cleaned_for_real_space_refine_molprobability.pdb

[command] rainbow

[menu] File -> Open -> .../user.map (set opacity into 50%)

[command] cd <user_path>/output/trajectory

[command] open traj.xtc structureModel #1

(note) Specified structure and coordinates should have the same number of atoms, if these are mismatched, I recommend to use Chimera instead

[command] coordset slider #1

Click black play button to play

Click red record button to record movie as .mp4 file

VMD (confirmed with version 1.9.3).

[menu] File -> New Molecule -> Browse -> (.../user.map) -> Load

[menu] Graphics -> Representation -> Draw -> Solid Surface

[menu] Graphics -> Representation -> Material -> Transparent

[menu] Graphics -> Representation -> Show -> Isosurface

[menu] File -> New Molecule -> Browse -> (output/trajectory/trajectory.gro) -> Load

Note

If the cryo_fit has been automatically repeated to improve cc, the trajectory file shows only the last cryo_fit run.

Therefore, to watch a full change of conformations, "phenix.cryo_fit user.mrc user.pdb no_rerun=True number_of_steps_for_cryo_fit=<enough step>" is recommended.