

## Install Cryo\_fit

### Contents

- [1. Install Phenix](#)
- [2. Link cryo\\_fit to phenix](#)
- [3. Install gromacs\\_cryo\\_fit](#)
- [4. Troubleshooting](#)

## 1. Install Phenix

See the [installation notes for Phenix](#)

## 2. Link cryo\_fit to phenix

At <user phenix>/modules

```
git clone git@github.com:cryoFIT/cryo_fit.git
```

At any folder,

```
libtbx.configure cryo_fit
```

## 3. Install gromacs\_cryo\_fit

Download gromacs-4.5.5\_cryo\_fit.zip from [here](#)

The screenshot shows the GitHub repository page for `cryoFIT / cryo_fit_install`. The repository has 2 watches, 0 stars, and 0 forks. The file `gromacs-4.5.5_cryo_fit.zip` is highlighted in the commit history, showing it was changed by `kimdn` 15 days ago. A blue arrow points from the file name to the `Download` button at the bottom right of the file view section.

Install using a provided script.

```
python <user phenix>/modules/cryo_fit/steps/0_install_cryo_fit/install_cryo_fit.py <cryo_fit.zip>
```

For example,

```
python /Users/doonam/bin/phenix-dev-2906/modules/cryo_fit/steps/0_install_cryo_fit/install_cryo_fit.py ~/Downloads/gromacs-4.5
```

## **4. Troubleshooting**

If user's OS doesn't have git

-> [install git](#)

If a user sees the error during git clone, "Permission denied (publickey).fatal: Could not read from remote repository. Please make sure you have the correct access rights and the repository exists.",

-> user's computer needs to be registered to his/her github account. Therefore, make github account (free), [Generate a new SSH key](#) and [Add a new SSH key to your GitHub account](#)

If user's macOS shows, "xcode-select: note: no developer tools were found at '/Applications/Xcode.app', requesting install. Choose an option in the dialog to download the command line developer tools."

```
download commandline as instructed (this may take ~7 minutes)
```

"-bash: libtbx.configure: command not found"

```
Run export PATH="/Users/.../phenix-xxxx/build/bin":$PATH
source /Users/.../phenix-xxxx/phenix_env.sh
in commandline, and try "libtbx.configure cryo_fit" again
```

"subprocess.CalledProcessError: Command '['which', 'phenix.cryo\_fit']'"

```
It can't find phenix.cryo_fit, link (sync) cryo_fit by "libtbx.configure cryo_fit".
Open a new tab in your terminal (like iterm2), then your .bash_profile should run freshly.
```