



## **Cryo\_fit1 FAQ**

### **Contents**

- [How can I provide gaussian filtered maps?](#)
- [How long does it take to run cryo\\_fit?](#)
- [How to enlarge map box size?](#)
- [How to extract a relevant map region?](#)
- [How to generate and record movie?](#)
- [How to improve initial cc?](#)
- [I can't run phenix.superpose\\_pdb with cryo\\_fitted\\_pdb file](#)
- [I observe cc values keep decreasing \(or do not increase that much\) during step 8](#)
- [I see "Fatal error: A charge group moved too far between two domain decomposition steps. This usually means that your system is not well equilibrated" at my 8\\_cryo\\_fit step.](#)
- [I see "Fatal error: Incomplete ring in HIS50" at my 8\\_cryo\\_fit step](#)
- [I see "Fatal error: Number of grid cells is zero. Probably the system and box collapsed." at my 8\\_cryo\\_fit step.](#)
- [I see "pdb file cleaning is not done" at my step 1 \(Make gro and topology file by regular gromacs\)](#)
- [I see "state.cpt not found, step 8 may be full of stepxb\\_nx.pdb" at my 8\\_cryo\\_fit step](#)
- [I see "step 0 correlation coefficient: nan" at my 8\\_cryo\\_fit step.](#)
- [I see "step 0 correlation coefficient: nan" and "Range checking error: Explanation: During neighbor searching, we assign each particle to a grid based on its coordinates. If your system contains collisions or parameter errors that give particles very high velocities you might end up with some coordinates being +-Infinity or NaN \(not-a-number\). Obviously, we cannot put these on a grid, so this is usually where we detect those errors. Make sure your system is properly energy-minimized and that the potential energy seems reasonable before trying again." at my 8\\_cryo\\_fit step.](#)
- [I see "The initial cell size \(xxx\) is smaller than the cell size limit \(xxx\), change options -dd, -rdd or -rcon, see the log file for details" at my 8\\_cryo\\_fit step](#)
- [I see "Too many LINCS warnings" at my 8\\_cryo\\_fit step](#)
- [I see "User's provided atomic model had 0.0 cc" in my cryo\\_fit.overall\\_log.](#)

### **How can I provide gaussian filtered maps?**

- See 'How can I provide gaussian filtered maps?' in [cryo\\_fit2\\_FAQ](#)


### **How long does it take to run cryo\_fit?**

# How long does it take to run cryo\_fit?

- CentOS with 2 cores (2.7 Ghz)
  - tRNA (6,000 atoms) : 2.5 minutes
  - Nucleosome (26,000 atoms): 1 hr
  - Beta-galactosidase (64,000 atoms): 2.5 hrs
  - ribosome (386,000 atoms): 3~7 hrs
- macOS
  - 24 cores (2.7 Ghz)
    - (4~16 cores would give similar/better performance)
    - Beta-galactosidase (64,000 atoms): 3.9 hrs
  - 4 cores (3.1 Ghz)
    - Beta-galactosidase (64,000 atoms, 50,000 cryo\_fit steps): 5.8 hrs



**Benchmark with ribosome** (same linux machine, 10k emsteps, 50k number\_of\_steps, 3x emweight)



- Based on this result, number of cores is recommended up to 16

## How to enlarge map box size?

- One can know map size by
  1. cryo\_fit which reports a map size between step 6 and step 7 with codes of
    - from iotbx import ccp4\_map
    - ccp4\_map = ccp4\_map.map\_reader(user\_input\_map)
    - target\_map\_data = ccp4\_map.map\_data()
  2. EMAN2's [e2iminfo](#)
    - For example, e2iminfo.py <user>.mrc
    - (note) As of EMAN2.2, the file extension should be mrc, not ccp4
  3. VMD which can visualize current map size as in this FAQ.
- Then, [relion\\_image\\_handler](#) can enlarge map box size.
  - For example, relion\_image\_handler --i <user>.mrc --new\_box 370 --o <user\_box\_size\_370>.mrc

## How to extract a relevant map region?

- Question: I only want to fit the monomer because only one monomer of the trimer map is the best. Do you have any idea of how to generate the monomer map that cryo\_fit can use?
- Answer: Please use [phenix.map\\_box](#) to extract local map.
  - Alternatively, you can try [phenix.map\\_to\\_model](#).

## **How to generate and record movie?**

- [Generate record movie by cryo\\_fit1](#)

## **How to improve initial cc?**

- [phenix.dock\\_in\\_map](#), UCSF Chimera's 'Fit in Map (Tools -> Volume Data -> Fit in Map)', UCSF Chimera's manual fitting (using a mouse) often improve initial cross-correlation (cc) between map and model. UCSF ChimeraX's isolate may improve initial cc as well.

## **I can't run phenix.superpose\_pdb with cryo\_fitted pdb file**

- phenix.superpose\_pdb can superimpose only between pdb files that have equally/similarly aligned nucleic names. This applies to all pdb input files (not only cryo\_fitted files). Consider to align both input files by "python <Phenix Path>/modules/cryo\_fit/steps/9\_after\_cryo\_fit/align\_nucleic\_acid\_name\_into\_middle/align\_nucleic\_acid\_name\_into

## **I observe cc values keep decreasing (or do not increase that much) during step 8**

- Initial few steps of cryo\_fit in step 8 may decrease cc values temporarily as an initial perturbation.
- However, if it keeps decreasing continuously (or do not increase that much), it may indicate two possible cases.
  - Case 1 <A user input pdb file has reasonable structural geometry> For this case, try below possible solutions. If cryo\_fit still can't find better (higher) cc, then the initial correlation between user input pdb file and cryo-EM map is already high enough. Just run phenix.real\_space\_refine only and deposit. Doo Nam recommends to argue/claim in a paper like "Our fit is so high, even cryo\_fit did not find higher level of fit than our initial fit".
  - Case 2 <A user input pdb file has unreasonable structural geometry> For this case, although the initial fit between user's atomic model and map looks good, it is a fictitious fitting without consideration of ideal atomistic model geometry. Run cryo\_fit to find decent fit to the cryo-EM map that restores/maintains reasonable structural geometry based on molecular dynamics forcefield (amber03.ff).
- There are possible solutions
  1. Provide a higher (better) resolution map tends to find higher cc.
  2. Cryo\_fit calculates the gradient of cc. If cryo\_fit is provided a giant cryo-em map with a tiny atomic model, then there is a large empty space (not filled). Therefore, the constraint forces for the model are very small. Consequently, these small forces are not helpful.
    - 2-1. Re-run cryo\_fit with an atomic model that fits the majority of the map. If you fit multiple atomic models into a symmetric map or do sequential fitting into a non-symmetric map, watch [Tom Goddard's lecture \(2013\)](#).
    - 2-2. Re-run cryo\_fit with only a relevant map region. See 'How to extract a relevant map region?' in this FAQ.
  3. If the initial model is not properly aligned to a map, see 'How to improve initial cc?' in this FAQ.
  4. A user may enforce a stronger initial map weight (e.g. emweight\_multiply\_by) manually. However, since cryo\_fit automatically increases emweight\_multiply\_by if it doesn't find higher cc anyway, manual increase in map weight may not necessarily find higher cc.

## **I see "Fatal error: A charge group moved too far between two domain decomposition steps. This usually means that your system is not well equilibrated" at my 8 cryo\_fit step.**

- Using macOS 10.13.6 helped rather than using Ubuntu 16.04. Maybe macOS has better numerical stability.
- One gromacs expert suggested to try smaller time\_step\_for\_cryo\_fit.
- However for Doo Nam, simply using macOS solved the problem.
- Less likely, but still a possible reason is that the map weight is too high. Therefore, lowering emweight\_multiply\_by may help.

## **I see "Fatal error: Incomplete ring in HIS50" at my 8 cryo\_fit step**

- Doo Nam observed this error when he tried to run cryo\_fit with 3jch.pdb that has incomplete Histidine ring atoms.
- He solved with UCSF Chimera

1. Open pdb file
  2. [menu]
  3. Select -> Residue -> HIS
  4. Tools -> Structure Editing -> Rotamers -> OK
  5. (select the most probable rotamer each)
  6. File -> Save PDB
- Pymol and Swissspdb viewer have similar residue fixing functions
  - Alternatively, a user may identify residues manually by running `<phenix-xxxx>/modules/cryo_fit/steps/0_prepare_cryo_fit/count_number_of_atoms_in_each_residue/`

### **I see "Fatal error: Number of grid cells is zero. Probably the system and box collapsed." at my 8 cryo fit step.**

- step\_8 may be full of stepxb\_nx.pdb.
  - Most likely, this means that initial cc is too low for MD simulation.
    - When Doo Nam ran `real_space_refine` first, then run `real_space_refined` atomic model in `cryo_fit`, it was solved.
    - Alternatively, improve initial cc by fitting initial atomic model into a map (see "How to improve initial cc?" in this FAQ)
  - Less likely, but still a possible case is when the map weight is too high, lowering `emweight_multiply_by` may help.

### **I see "pdb file cleaning is not done" at my step 1 (Make gro and topology file by regular gromacs)**

"I edited out lipids, HEM and other hetero atoms and I verified that they are all gone. However, still my pdb file is not clean enough for `gromacs_cryo_fit`".

- If your lipids, HEM and other hetero atoms are not necessary for now, you may take out as you did. You may add those atoms after flexible fit (`cryo_fit`). Please send Doo Nam ([doonam@lanl.gov](mailto:doonam@lanl.gov)) your pdb input file only (no map). He will take a look. `gromacs_cryo_fit` uses `amber03` which may not have forcefield parameters for your less-usual atoms/residues.
- If your lipids, HEM and other hetero atoms are necessary for now and your target is protein, Doo Nam would recommend [cryo\\_fit2](#) since it supports `phenix.eLBOW` derived ligand parameters.

### **I see "state.cpt not found, step 8 may be full of stepxb nx.pdb" at my 8 cryo fit step**

- Literally, `steps/8_cryo_fit` fold may be full of `stepxb_nx.pdb` files.
  - Most likely, this means that a `cryo_fit` input pdb file is not yet stable enough for sensitive Gromacs MD simulation (step 8) even after `cryo_fit`'s ample minimization step (e.g. step 4).
    - When Doo Nam ran `real_space_refine` first, then run `real_space_refined` atomic model in `cryo_fit`, it was solved more than 2 cases.
    - If initial cc is too low, improve initial cc by fitting initial atomic model into a map,
      - either by following "How to improve initial cc?" in this FAQ
      - or adding/copying more atoms.
        - For example, when Doo Nam tried to fit a monomer into a trimer map, this error occurred. However, simply adding dimer atomic models into an input pdb file and running `cryo_fit` again removed the error.
    - Less likely, but still a possible case is when the map weight is too high, lowering `emweight_multiply_by` may help.


### **I see "step 0 correlation coefficient: nan" at my 8 cryo fit step.**

- This often indicates that the initial atomic model is not placed into a cryo-EM map.
  - Therefore, please improve initial cc by fitting initial atomic model into a map (see "How to improve initial cc?" in this FAQ)
  - If a user observed "Range checking error: Explanation:" error message at 8\_cryo\_fit step, please refer `<I see "step 0 correlation coefficient: nan" and "Range checking error: Explanation:">` part in this FAQ


### **I see "step 0 correlation coefficient: nan" and "Range checking error: Explanation: During neighbor searching, we assign each**

**particle to a grid based on its coordinates. If your system contains collisions or parameter errors that give particles very high velocities you might end up with some coordinates being +- Infinity or NaN (not-a-number). Obviously, we cannot put these on a grid, so this is usually where we detect those errors. Make sure your system is properly energy-minimized and that the potential energy seems reasonable before trying again." at my 8 cryo fit step.**

- Based on Doo Nam's experiences, this error message may indicate between two cases/scenarios.
  1. Doo Nam observed this error when starting molecule is not stabilized structurally. There are two possible solutions.
    - Run phenix.real\_space\_refine before cryo\_fit. Then, provide real\_space\_refined molecule into cryo\_fit. Then, most problems were solved.
    - When Doo Nam skipped minimization, this error occurred. Since most users minimize the starting structure by default (cryo\_fit default), most users do not need care about this. Besides, as of 9/18/2019, the default value of number\_of\_steps\_for\_minimization=20000 (which is a lot).
  2. It means that the map dimensions need to be larger. Therefore, check map size and solve a problem with VMD/EMAN2/relion like followings.
    - Like other MD simulations, gromacs need enough map box size to cover the atomic model to run (ziggle and wiggle). Refer [Waters seems to be out of the box](#)
    - For example, stuck-out red oxygen atoms outside the right edge of the box are the problem.



- In order to run any MD simulation (including cryo\_fit), a box should be large enough like



- Make map box size larger (see "How to enlarge map box size?" in this FAQ), and run cryo\_fit again. You can check map box size by VMD. Alternatively, remove sticking out atoms if these are unnecessary, then run cryo\_fit again.
- For protein modeling, I would use [cryo\\_fit2](#) which is not limited by box size requirement. Most of the time, it better fits than cryo\_fit1 in terms of fitting and geometry preservation anyway.

**I see "The initial cell size (xxx) is smaller than the cell size limit (xxx), change options -dd, -rdd or -rcon, see the log file for details" at my 8 cryo\_fit step**

- Doo Nam observed this error for 2 cases
  1. When he provided a pdb file that has unspecified/deleted region of a molecule. Even when vmd assisted map box/cell dimension is larger than initial molecule space, this error

appeared.

2. When he provided a pdb file that has small number of residues (e.g. 5~16 amino acids). Even when vmd assisted map box/cell dimension is larger than initial molecule space, this error appeared. This error appeared regardless of the existence of CRYST1 and SCALE header.

### **I see "Too many LINCS warnings" at my 8 cryo fit step**

- Doo Nam observed this error when starting molecule is not stabilized structurally. There are 3 possible solutions.
  1. Run phenix.real\_space\_refine before cryo\_fit. Then, provide real\_space\_refined molecule into cryo\_fit. Then, most problems were solved.
  2. When Doo Nam skipped minimization, this error occurred. Since most users minimize the starting structure by default (cryo\_fit default), most users do not need care about this. Besides, as of 9/18/2019, the default value of number\_of\_steps\_for\_minimization=20000 (which is a lot).
  3. If starting molecule has non-standard chemical geometry that is too extreme that are not stabilized neither by real\_space\_refine nor by cryo\_fit minimization, then a user needs to fix those non-standard geometry manually.

### **I see "User's provided atomic model had 0.0 cc" in my cryo\_fit.overall log.**

- It means that literally user's provided atomic model has no correlation to map to start with.
- See "How to improve initial cc?" in this FAQ, and provide initially aligned atomic model to cryo\_fit.